



Tech & Firmware Notification  
PXE/USB Imaging Utility 5.22.0110

January 31, 2021



---

10ZiG are pleased to announce a new version of their PXE and USB Imaging Utility, now featuring support for 10ZiG Linux NOS™ 64 Bit based compressed backup and restore. Due to increasing storage capacity and image sizes now available, the compressed method offers improved scalability & compatibility for performing backup and restore of 10ZiG Linux NOS™ images across applicable 10ZiG hardware.

Customer's and Partner's are highly encouraged to upgrade to this version when performing backup or restore of 10ZiG Linux NOS™ images via the 10ZiG Manager (PXE) or USB Imaging Utility.

Refer to the 10ZiG imaging Utility Guide for further guidance on usage.

### **Compatible 10ZiG Hardware**

4400 - 4500 - 4600q - 5800q - 5800qd - 5900q - 6000q - 6100

### **Compatible 10ZiG Firmware**

NOS™ - all versions

### **Feature Requests Implemented**

10ZiG Linux NOS™ Backup process switched from full disk based binary (BIN) based to Compressed method (BIM).

10ZiG Linux NOS™ Restore process compatible with Compressed method (TAR) and partitions storage accordingly.

### **Bug Fixes**

Cannot mount network share with spaces in the path.

Intermittently connecting to a network image (SMB) fails most commonly on the 1st attempt.